



Features and Benefits

- **EM sensor technology**
<1e⁻ read noise
- **Fringe suppression technology as standard (970-BVF only)**
Fringing minimized for NIR applications
- **Multi-Megahertz Readout**
High repetition rates achievable with low noise electronics
- **TE cooling to -100°C**
Negligible dark current without the inconvenience of LN₂
- **UltraVac™**
Critical for sustained vacuum integrity and to maintain unequalled cooling and QE performance, year after year
- **16 x 16 μm pixel size**
Optimized pixel size for high resolution spectroscopy
- **Dual output amplifiers**
Software-selectable between conventional High Sensitivity output (low light) or an Electron Multiplying output (ultra low light)
- **Crop Mode operation**
Achieve the highest possible spectral rates of over 1,500 spectra per second
- **USB 2.0 connection**
Ideal for laptop operation
Seamless operation alongside USB-based Shamrock spectrograph family
- **Solis software for Spectroscopy**
Comprehensive, user-friendly interface for simultaneous detector & spectrograph control
- **Software Development Kit (SDK)**
Ease of control integration into complex setups: Matlab, Labview, Visual Basic or C/C++

Market leading platform for ultrasensitive and ultrafast spectroscopy.

EM technology enables charge from each pixel to be multiplied on the sensor before readout, providing single photon sensitivity. The Newton EM platform combines a 1600 x 200 (or 1600 x 400) array of 16 μm pixels, thermoelectric cooling down to -100°C for negligible dark current, 3MHz readout and USB 2.0 plug-and-play connectivity to provide unrivalled performance for spectroscopic applications. The dual output amplifiers allow software selection between either a conventional High Sensitivity or Electron Multiplying outputs to suit a broad range of photon regime conditions. This makes the Newton EMCCD the ideal choice for ultrafast chemical mapping applications e.g. SERS, TERS or luminescence mapping.

Specifications Summary ^{*1}

Active pixels	1600 x 200 or 1600 x 400
Pixel size (W x H)	16 x 16 μm
Image area	25.6 x 3.2 or 6.4 mm
Output node well depth (typical)	
High Sensitivity mode	300,000 e ⁻
Electron Multiplying mode	1,300,000 e ⁻
Maximum cooling	-100°C
Maximum spectra per sec	1,515
Read noise	As low as 2.8 e ⁻ (< 1 e ⁻ in EM mode)
Dark current	As low as 0.00007 e ⁻ /pixel/sec

Key Specifications ^{•1}

Model number	DU970P	DU971P
Sensor options	• BV: Back Illuminated CCD, Vis-optimized • BVF: Back Illuminated CCD, Vis-optimized and anti-fringing • FI: Front Illuminated CCD • UV: Front Illuminated CCD with UV coating • UVB: Back Illuminated CCD with UV coating	• BV: Back Illuminated CCD, Vis-optimized • FI: Front Illuminated CCD • UV: Front Illuminated CCD with UV coating • UVB: Back Illuminated CCD with UV coating
Active pixels ^{•2}	1600 x 200	1600 x 400
Pixel size	16 x 16 μm	
Image area	25.6 x 3.2 mm with 100% fill factor	25.6 x 6.4 mm with 100% fill factor
Minimum temperatures ^{•3}		
Air cooled	-80°C	
Coolant recirculator	-95°C	
Coolant chiller, coolant @ 10°C, 0.75l/min	-100°C	
Max spectra per second ^{•4}	649 (Full Vertical Bin), 1,515 (Crop Mode - 20 rows)	396 (Full Vertical Bin), 1,515 (Crop Mode - 20 rows)
System window type	Single UV-grade fused silica window, uncoated. Various AR coatings & MgF ₂ options available	
Blemish specifications	Grade 1 as per sensor manufacturer definition	

Advanced Specifications ^{•1}

Dark current, e ⁻ /pixel/sec @ max cooling						
FI, UV	0.00007					
BV, UVB	0.00020					
BVF	0.00010					
Output node well depth (typical)						
High Sensitivity mode	300,000 e ⁻					
Electron Multiplying mode	1,300,000 e ⁻					
Register well depth						
High Sensitivity mode	400,000 e ⁻					
Electron Multiplying mode	800,000 e ⁻					
Active area pixel well depth	210,000 e ⁻ ^{•5}					
Read noise (e ⁻) ^{•6}	50 kHz	1 MHz	3 MHz	50 kHz	1 MHz	3 MHz
High Sensitivity mode: Typ (Max) - EM off	2.8 (5)	6.7 (9)	8.5 (12)	2.8 (5)	6.7 (9)	8.5 (12)
Electron Multiplying mode: Typ (Max) - EM off	8 (15)	25 (35)	38 (50)	8 (15)	25 (35)	38 (50)
Electron Multiplying mode: Typ (Max) - EM on	< 1	< 1	< 1	< 1	< 1	< 1
Sensitivity (e ⁻ /count)						
High Sensitivity mode	Adjustable from 0.8 - 3					
Electron Multiplying mode	Adjustable from 5 - 20					
Electron Multiplier gain	1 - 1,000 times (software controlled)					
Linearity ^{•6}	Better than 99%					
Digitization	16 bit					
Vertical clock speed ^{•8}	4.9, 9.8, 19, 38, 57 (software selectable)					


Have you found what you are looking for?

Need to work further into the NIR? The iDus InGaAs series, with up to 1024 pixel linear array with transmission to 2.2 μm .

Need a customized version? Please contact us to discuss our Customer Special Request options.


The Newton series combines seamlessly with Andor's research grade Shamrock Czerny-Turner spectrographs.

Dark Current ^{*9}




Quantum Efficiency Curves ^{*10}


25°C



Typical Setup



Readout Rate & Speed ^{*11}



Creating The Optimum Product for You

How to customize the Newton EMCCD series:

Step 1.

Simply select from the 2 sensor array types that best suit your needs from the selection opposite.

Step 2.

The Newton EMCCD comes with 5 options for sensor types. Please select the sensor which best suits your needs.

Step 3.

Please indicate if you wish to select an alternative window and which software you require.

Step 4.

For compatibility, please indicate which accessories are required.

DU 970 P- BVF

example shown

Step 1.

Choose sensor array

970P: 1600 x 200 pixel array
971P: 1600 x 400 pixel array

Step 2.

Choose sensor type

BV: Back Illuminated CCD, Vis-optimized
BVF: Back Illuminated CCD, Vis-optimized and anti-fringing (970 model only)
FI: Front Illuminated CCD
UV: Front Illuminated CCD with UV coating
UVB: Back Illuminated CCD with UV coating

Step 3.

The Newton EMCCD series models are supplied with an uncoated UV-grade fused silica window as standard. The following alternative window choices are available and must be ordered at time of build (if selected):

OPTION-C1-AR1 AR coated UV-grade fused silica window (optimized broadband visible 400-900nm). 50% transmission at 180 nm

OPTION-C1-MGF2 Magnesium Fluoride window for transmission in the VUV.
50% transmission at 120 nm

The Newton EMCCD also requires at least one of the following software options:

Solis for Spectroscopy A 32-bit application compatible with 32 and 64-bit Windows (XP, Vista, 7 and 8) offering rich functionality for data acquisition and processing. AndorBasic provides macro language control of data acquisition, processing, display and export. Control of Andor Shamrock spectrographs and a very wide range of 3rd party spectrographs is also available, see list below.

Andor SDK A software development kit that allows you to control the Andor range of cameras from your own application. Available as 32 and 64-bit libraries for Windows (XP, Vista, 7 and 8) and Linux. Compatible with C/C++, C#, Delphi, VB6, VB.NET, LabVIEW and Matlab.

Step 4.

The following accessories are available:

XW-RECR Coolant re-circulator for enhanced cooling performance

ACC-XW-CHIL-160 Oasis 160 Ultra Compact Chiller Unit (tubing to be ordered separately)

ACC-6MM-TUBING-2xxxM 6 mm tubing option for ACC-XW-CHIL-160

LM-C C-mount lens adaptor

LM-NIKON-F C-mount lens adaptor

LMS-NIKON-F-NS25B Nikon F-mount lens adaptor with shutter

ACC-SD-VDM1000 Shutter Driver for NS25B Bistable Shutter (not needed for Shamrock spectrographs)

ACC-SHT-NS25B Bistable Shutter, Standalone (not needed for Shamrock spectrographs)

Spectrograph Compatibility




The Newton series is fully compatible with Andor's Shamrock spectrograph (163 - 750 nm focal lengths) family. Spectrograph mounting flanges and software control are available for a wide variety of 3rd party spectrographs including, McPherson, JY/Horiba, PI/Acton, Chromex/Bruker, Oriel/Newport, Photon Design, Dongwoo, Bentham, Solar TII and others.



Newton EMCCD mounted on a Shamrock 750 mm triple grating imaging spectrograph, ideal for high resolution spectroscopy.

Product Drawings

Dimensions in mm [inches]



■ = position of pixel 1,1

Weight: 2.7 kg [5 lb15 oz]

Connecting to the Newton

Camera Control

Connector type: USB 2.0

TTL / Logic

Connector type: SMB, provided with SMB - BNC cable

1 = Fire (Output), 2 = External Trigger (Input), 3 = Shutter (Output)


I²C connector

Compatible with Fischer SC102A054-130


1 = Shutter (TTL), 2 = I²C Clock, 3 = I²C Data, 4 = +5 Vdc, 5 = Ground

Minimum cable clearance required at rear of camera

90 mm



Mounting hole locations



Rear connector panel

Applications & Techniques Guide

	BV models	BVF models	FI models	UV models	UVB models
Absorption/Transmittance/Reflection	✓	✓		✓	✓
Atomic Emission Spectroscopy	✓	✓	✓	✓	
Fluorescence & Luminescence	✓	✓	✓	✓	✓
Raman Spectroscopy (244 – 488 nm)	✓	✓		✓	✓
Raman Spectroscopy (514, 532 nm)	✓	✓		✓	
Raman Spectroscopy (633 nm)	✓	✓	✓		
Photon Counting	✓	✓			✓
Single Molecule Spectroscopy	✓	✓		✓	✓

✓ = Suitable

✓ = Optimum



Order Today

Need more information? At Andor we are committed to finding the correct solution for you. With a dedicated team of technical advisors, we are able to offer you one-to-one guidance and technical support on all Andor products. For a full listing of our local sales offices, please see: andor.com/contact

Our regional headquarters are:

Europe

Belfast, Northern Ireland
Phone +44 (28) 9023 7126
Fax +44 (28) 9031 0792

Japan

Tokyo
Phone +81 (3) 3518 6488
Fax +81 (3) 3518 6489

North America

Connecticut, USA
Phone +1 (860) 290 9211
Fax +1 (860) 290 9566

China

Beijing
Phone +86 (10) 5129 4977
Fax +86 (10) 6445 5401

Items shipped with your camera:

- 1x 2m BNC - SMB connection cable
- 1x 3m USB 2.0 cable Type A to Type B
- 1x Set of Allen keys (7/64", 3/32" & 3 mm)
- 1x Power supply with mains cable
- 1x Quick launch guide
- 1x CD containing Andor user guides
- 1x Individual system performance booklet
- 1x CD containing either Solis software or SDK (if ordered)

Footnotes: Specifications are subject to change without notice

1. Figures are typical unless otherwise stated.
2. Edge pixels may exhibit a partial response.
3. Cooling is provided by the use of an external mains driven power supply. Minimum temperatures listed are typical values with ambient temperature of 20°C. Systems are specified in terms of minimum dark current achievable rather than absolute temperature.
4. Based on horizontal pixel readout rate of 3 MHz and a vertical shift speed (in conventional mode) of 4.9 μ s. Achievable spectral rates will vary with selected trigger mode.
5. Shown for EM mode. For Conventional mode the measurable well depth value will be lower, as a result of the combination of higher sensitivity values and A/D 16 bits digitization.
6. Readout noise is for the entire system. It is a combination of CCD readout noise and A/D noise. Measurement is for Single Pixel readout with the CCD at a temperature of -80°C and minimum exposure time under dark conditions. Noise values will change with readout mode.
7. Linearity is measured from a plot of counts vs exposure time under constant photon flux up to the saturation point of the system.
8. Vertical speeds are software selectable. All sensors are designed to give optimum Charge Transfer Efficiency (CTE) at 9.7 μ s vertical pixel shift, some decrease in CTE may be observed at faster shift speeds.
9. The graph shows typical dark current level as a function of temperature. The dark current measurement is averaged over the CCD area excluding any regions of blemishes.
10. Quantum efficiency of the sensor at 20°C as supplied by the sensor manufacturer.
11. The chart shows the maximum possible readout rates available when using Multi-track mode, each track being defined as 20 rows. Crop mode is a specific single-track readout method optimized for rapid kinetic-type acquisition.

Minimum Computer Requirements:

- 3.0 GHz single core or 2.4 GHz multi core processor
- 2 GB RAM
- 100 MB free hard disc to install software (at least 1 GB recommended for data spooling)
- USB 2.0 High Speed Host Controller capable of sustained rate of 40 MB/s
- Windows (XP, Vista, 7 and 8) or Linux

Operating & Storage Conditions

- Operating Temperature: 0°C to 30°C ambient
- Relative Humidity: <70% (non-condensing)
- Storage Temperature: -25°C to 50°C

Power Requirements

- 110 - 240 VAC, 50 - 60 Hz



Windows is a registered trademark of Microsoft Corporation.
Labview is a registered trademark of National Instruments.
Matlab is a registered trademark of The MathWorks Inc.

SNewtonEMSS 0713 R1