

### 3.7 What is $M^2$ ?

$M^2$  or Beam Propagation Ratio, is a value that indicates how close a laser is to being a single mode  $TEM_{00}$  beam, which in turn determines how small a beam waist can be focused. For the perfect Gaussian  $TEM_{00}$  condition the  $M^2$  equals 1.

For a laser beam propagating through space, the equation for the divergence  $\theta_0$  of an unfocused beam is given by:

$$\theta_0 = M^2 4\lambda / \pi D_0 \quad (D_0 \text{ is the waist diameter of the laser beam})$$

For a pure Gaussian  $TEM_{00}$  beam  $M^2$  equals 1, and thus has no impact on the calculation. The calculation of the minimal beam spot after the lens is then:

$$d_0 = 4\lambda / \pi \theta \quad (\theta \text{ is the beam divergence after the lens})$$

Again with  $M^2$  equal to 1, the focused spot is diffraction limited. For real beams,  $M^2$  will be greater than 1, and thus the minimum beam waist will be larger by the  $M^2$  factor.




Figure 1 - Characteristics of a laser beam as it passes through a focusing lens

### How is $M^2$ measured?

$M^2$  cannot be determined from a single beam profile measurement. The ISO/DIS 11146 requires that  $M^2$  be calculated from a series of measurements as shown in figure 1.  $M^2$  is measured on real beams by focusing the beam with a fixed position lens of known focal length, and then measuring the characteristics of the artificially created beam waist and divergence.

To provide an accurate calculation of  $M^2$ , it is essential to make at least 5 measurements in the focused beam waist region, and at least 5 measurements in the far field, two Rayleigh ranges away from the waist area. The multiple measurements ensure that the minimum beam width is found. In addition, the multiple measurements enable a “curve fit” that improves the accuracy of the calculation by minimizing measurement error at any single point. An accurate calculation of  $M^2$  is made by using the data from the multiple beam width measurements at known distances from a lens, coupled with the known characteristics of the focusing lens.

### $M^2$ Measurement Solutions

Ophir-Spiricon have a number of solutions for the measurement of  $M^2$  ranging from simple manual processes to fully automated dedicated instruments, depending on the frequency of the need to measure  $M^2$  of lasers and laser systems. We have a system that will meet most needs, whether for research and development of new laser systems, manufacturing quality assurance, or maintenance and service of existing systems.




Figure 2 - focusing laser beam